

## **Optimization of Real Time Quantitative PCR (Q-PCR) for *Fusarium pseudograminearum* and *F. culmorum* on wheat**

G. J. POOLE (1), F. Ozdemir (2), S. D. Nydam (1), K. L. Schroeder (3), T. C. Paulitz (3), J. M. Nicol (4), K. G. Campbell (5)

(1) Washington State University, Pullman, WA, USA; (2) Bahri Dagdas International Agricultural Research Institute, Konya, Turkey; (3) USDA-ARS Root Disease and Biological Control Research Unit; (4) CIMMYT (International Maize and Wheat Improvement Centre) Emek 06511 Ankara, Turkey; (5) USDA-ARS Wheat Genetics, Quality, Physiology and Disease Research Unit

Phytopathology 99:S103

*Fusarium* crown rot of wheat is caused by a complex of *Fusarium* species, of which *F. pseudograminearum* and *F. culmorum* are the most important. Crown rot reduces wheat yields by an average of 9% in the Pacific Northwest. Traditional methods of species identification have included morphological characteristics of macroconidia. With the advent of Q-PCR techniques and the development of primers for *F. pseudograminearum* and *F. culmorum*, the potential exists for more accurate species identification and fungal DNA quantification from infected wheat stems. Primers developed in previous studies were evaluated for use in Q-PCR and DNA extraction kits were tested and optimized to accurately assess *Fusarium* species and DNA concentrations in wheat stem tissue. The 'OPT' primers (Shilling et al. 1996) for the amplification of *F. culmorum* and the 'FPG' primers (Williams et al. 2002) for the amplification of *F. pseudograminearum* yielded the most consistent results. The MO-BIO® Ultra Clean Soil Kit for DNA extraction produced the most consistent Q-PCR amplification from infected wheat tissue. The most optimal results were obtained by grinding with liquid nitrogen and soaking prior to bead beating (using a ceramic bead (MP Biomedicals® - FastDNA extraction kit)) and a Fast Prep speed of 5 for 45s. Addition of polyvinylpyrrolidone (PVPP) was necessary for adequate DNA extraction.